

FACE MASKS

MN 126 | Disposable Anesthesia Face Mask

Tek Kullanımlık Anestezi Maskesi

DISPOSABLE FACE MASK

CODE	SIZE	DEAD SPACE
MN 126-01	00 Pre-neonate	20 ml.
MN 126-02	0 Neonate	25 ml.
MN 126-03	1 Infant	45 ml.
MN 126-04	2 Child	85 ml.
MN 126-05	3 Small	180 ml.
MN 126-06	4 Medium	220 ml.
MN 126-07	5 Large	250 ml.
MN 126-08	6 Ultra Large	310 ml.



MN 127 | Reusable Anesthesia Face Mask

Çok Kullanımlık Anestezi Maskesi



Black Rubber

REUSABLE BLACK RUBBER

CODE	SIZE
MN 127-01	000
MN 127-02	00
MN 127-03	0
MN 127-04	1
MN 127-05	2
MN 127-06	3
MN 127-07	4
MN 127-08	5
MN 127-09	6



Black Rubber



Silicone

REUSABLE SILICONE

CODE	SIZE
MN 127-10	0
MN 127-11	1
MN 127-12	2
MN 127-13	3
MN 127-14	4
MN 127-15	5



Semi-Silicone

REUSABLE SEMI - SILICONE

CODE	SIZE
MN 127-16	3
MN 127-17	4
MN 127-18	5

Reusable face masks can be sterilized (autoclavable) up to 134 °C without any deformation.

Çok kullanımlık anestezi maskeleri 134 °C kadar otoklavda herhangi bir bozulma olmadan steril edilebilir.