

MN 131 Reusable Manual Artificial Resuscitation Bags

Çok Kullanımlık Suni Solunum Balonları



Manual resuscitator (Silicone). (including silicone bag, PSF reservoir valve set, PSF intake valve set, silicone mask 2 mt. oxygen tubing, PVC reservoir bag, 60 cm H₂O pop off.), Mouth opener, 3 different size airway

CODE	SIZE
MN 131-01	NEONATE
MN 131-02	INFANT
MN 131-03	PEDIATRIC
MN 131-04	ADULT



SILICONE



Airways

Mouth opener

Reusable Silicone and Black Rubber manual resuscitators can be sterilized (autoclavable) up to 134 °C without any deformation.
Çok kullanımlık Silikon ve kauçuk suni solunum balonları 134 °C kadar otoklavda herhangi bir bozulma olmadan steril edilebilir.

MN 131 Disposable Manual Artificial Resuscitation Bags

Tek Kullanımlık Suni Solunum Balonları

Manual resuscitator (PVC). (including silicone bag, PSF reservoir valve set, PSF intake valve set, PVC mask 2 mt. oxygen tubing, PVC reservoir bag, 60 cm H₂O pop off.)

CODE	SIZE
MN 131-05	ADULT
MN 131-06	PEDIATRIC



PVC