

MN 138 | Mortonvent-T Tube and Filter
Mortonvent-T Parçası ve Filtresi



MN 138-01



MN 138-03

CODE	MODEL
MN 138-01	MORTONVENT SET (Filter&Oxygen Tube)
MN 138-03	MORTONVENT FILTER WITH OXYGEN PORT

TECHNICAL PROPERTIES

MORTONVENT and accessory products, the MORTONVENT Holder is designed specifically for trachestomized or intubated spontaneously breathing patients.

The MORTONVENT weight is only 4.0 grams and has a dead space of just 10 ml. Recommended for tidal volumes of 50 ml. and above MORTONVENT can be used on both pediatric and adult patients.

Mortonvent holder conveniently attaches MORTONVENT filter to tracheostomy or endotracheal tube, allowing it to easily be removed and replaced during suctioning or change out MORTONVENT filter to allow delivery of supplemental oxygen.