

MN 138 Mortonvent-T Tube and Filter

Mortonvent-T Parçası ve Filtresi



MN 138-01



MN 138-03

CODE	MODEL
MN 138-01	MORTONVENT SET (Filter&Oxygen Tube)
MN 138-03	MORTONVENT FILTER WITH OXYGEN PORT

TECHNICAL PROPERTIES

MORTONVENT and accessory products, the MORTONVENT Holder is designed specifically for tracheostomized or intubated spontaneously breathing patients.


The MORTONVENT weight is only 4.0 grams and has a dead space of just 10 ml. Recommended for tidal volumes of 50 ml. and above MORTONVENT can be used on both pediatric and adult patients.

Mortonvent holder conveniently attaches MORTONVENT filter to tracheostomy or endotracheal tube, allowing it to easily be removed and replaced during suctioning or change out MORTONVENT filter to allow delivery of supplemental oxygen.

MN 139 Disposable Laryngeal Airways

Tek Kullanımlık Laryngeal Airways


Eksalasyon Valfli, Tek Su Tutuculu Solunum Devresi



CODE	SIZE	PATIENT (kg)	MAX. CUFF AIR VOL.
MN 139-01	DISPOSABLE SIZE:1	5 kg.	4.0 ml.
MN 139-02	DISPOSABLE SIZE:1,5	5 - 10 kg.	7.0 ml
MN 139-03	DISPOSABLE SIZE:2	10 - 20 kg.	8.0 ml
MN 139-04	DISPOSABLE SIZE:2,5	20 - 30 kg.	11.5 ml
MN 139-05	DISPOSABLE SIZE:3	30 - 50 kg.	20 ml
MN 139-06	DISPOSABLE SIZE:4	50 - 70 kg.	30 ml
MN 139-07	DISPOSABLE SIZE:5	70 - 100 kg.	30 ml

MN 140 Reusable Laryngeal Airways

Çok Kullanımlık Laryngeal Airways



CODE	SIZE	PATIENT (kg)	MAX. CUFF AIR VOL.
MN 140-01	REUSABLE SIZE:1	5 kg.	4.0 ml.
MN 140-02	REUSABLE SIZE:1,5	5 - 10 kg.	7.0 ml
MN 140-03	REUSABLE SIZE:2	10 - 20 kg.	8.0 ml
MN 140-04	REUSABLE SIZE:2,5	20 - 30 kg.	11.5 ml
MN 140-05	REUSABLE SIZE:3	30 - 50 kg.	20 ml
MN 140-06	REUSABLE SIZE:4	50 - 70 kg.	30 ml
MN 140-07	REUSABLE SIZE:5	70 - 100 kg.	30 ml