

**MN 143** | Nebulizer Bottle  
Nebulizatör Haznesi



**TECHNICAL PROPERTIES**

NEBULIZER BOTTLE
Material: PP (cap), PS (inner cone, bottle)
Option color available, min. order quantities request.
Feature: 18mm OD connector,
Nebulizes 3cc within 7 minutes 8LPM 6cc volume

**MN 143 - 01** | Valved Tee Adapter  
Valfli T Konnektör



**TECHNICAL PROPERTIES**

15M / 15F
-----------

**MN 143-02** | Inhalation Set with T Valved  
T Valfli İnhalasyon Seti



**TECHNICAL PROPERTIES**

JET NEBULIZER
T PIECE 22 M / 15 F
OXYGEN TUBE
CORRUGATED TUBE

**MN 143-03** | Inhalation Chamber with T Valved  
T Valfli İnhalasyon Tedavi Çemberi



**TECHNICAL PROPERTIES**

JET NEBULIZER
T PIECE 10 M / 15 F
OXYGEN TUBE

**MN 144** | Jet Set Nebulizer  
Jet Nebulizatör Seti



**TECHNICAL PROPERTIES**

JET NEBULIZER
T PIECE 22 M / 15 M
OXYGEN TUBE

**MN 144 - 01** | Jet Set Nebulizer with mouth piece  
Ağızlıklı Jet Nebulizatör Seti



**TECHNICAL PROPERTIES**

JET NEBULIZER
T PIECE
MOUTH PIECE
OXYGEN TUBE