

MN 147 | **Triball-Breathing Exercises System**
Üç Top -Solunum Egzersiz Cihazı



TECHNICAL PROPERTIES

% 97 CALIBRATION

With the specially developed intermediary connector part measurement is made with 97% precision. Each product is tested for precision calibration.

The device consists of three graduated cylinders for indicating flow rate.

Different coloured balls in the cylinders indicate the approximate flow passing through the device.

A tube with mouth-piece is also included.

TEKNİK ÖZELLİKLER

Geliştirilen özel ara konektör ile %97 doğru ölçüm sağlanmıştır.

Her ürüne hassas kalibre testi yapılmaktadır.