

MN 135 | **Humidifier Chamber**
Nemlendirici Hazne

Adult



TECHNICAL PROPERTIES

Compressible volume full 260 ml,
Empty capacity 560 ml, max.
Water capacity 300 ml,
Reusable and disposable models are available
Inlet and outlet gas port 22 mm
Max. peak flow 180 L/min.
Materials used in construction :
Aluminium (Covered electroplate base), Polistiren (Dome), PVC Gasket)
It is a single use UV irradiated humidification chamber.
It is used to provide humidifaction for intubated or non-intubated patients.

CODE	SIZE
MN 135-01	HUMIDIFIER CHAMBER ADULT
MN 135-02	HUMIDIFIER CHAMBER ADULT AUTOFEED