

**MORTON**

**MN 400 | Pleural Drainage System 2000 cc.**  
Göğüs Drenaj Seti 2000 cc.



**TECHNICAL PROPERTIES**

2000 ml. collecting capacity.

Special design (prevent leaking).

Firm stand to avoid accidental overturn.

Easy write on drainage marking scale.

Handle for patient.

Length of tubing 150 cm.

Soft and kink resistance silk finished tubing.

Tube connector specially designed to adaptation different thoracic catheter sizes.

Luer lock port cap on the connector helps to empty coagulum collected in the tubing.

Sterile.

**TEKNİK ÖZELLİKLER**

2000 ml. kapasiteli.

Sızma riskini yok eden tek parça özel tasarım.

Devrilmeyi önleyen geniş taban.

Seviye göstergesi.

Pratik taşıma askısı.

Hortum uzunluğu 150 cm.

Bükülmeye dirençli, kırılmayan özel hortum.

Spiral hortum destekli.

Değişik boyutlardaki toraks kataterlerine bağlantı sağlayan portlu konnektör.

Konnektör üzerinde hava boşaltma kapağı.

Steril.