

MORTON

MN 401 | Pleural Drainage System 700 Cc.
Göğüs Drenaj Seti 700 Cc.



TECHNICAL PROPERTIES

- 700 ml. collecting capacity.
- Special design (prevent leaking).
- Firm stand to avoid accidental overturn.
- Easy write on drainage marking scale.
- Handle for patient.
- Length of tubing 150 cm.
- Soft and kink resistance silk finished tubing.
- Tube connector specially designed to adaptation different thoracic catheter sizes.
- Luer lock port cap on the connector helps to empty coagulum collected in the tubing.
- Sterile.

TEKNİK ÖZELLİKLER

- 700 ml. kapasiteli.
- Sızma riskini yok eden tek parça özel tasarım.
- Devrilmeyi önleyen geniş taban.
- Seviye göstergesi.
- Pratik taşıma askısı.
- Hortum uzunluğu 150 cm.
- Bükülmeye dirençli, kırılmayan özel hortum.
- Spiral hortum destekli.
- Değişik boyutlardaki toraks kataterlerine bağlantı sağlayan portlu konnektör.
- Konnektör üzerinde hava boşaltma kapağı.
- Steril.