

**MN 205 | I.V Regulator**  
Damla Ayar Seti



**TECHNICAL PROPERTIES**

**Tolerance of flow rate;**

5-20 ml/h :30%

20-40 ml/h :20%

40-250 ml/h :15%

**Operating range:**

5 to 250 ml/h

**Diameter x Height:**

32.00 x 34.20mm. (body)

Inlet/ Outlet connectors

Weight:9.00 gr.

Standard IV tubing connectors (2.7x 3.5 or 3.0 x 4.1 mm)



***New Model***